



Cushioned top edge for patient comfort

Patented torso wrap design provides patient stability and reduces risk of skin tears, bruising and fractures.

19-25 handles allow for multiple caregivers to team lift while using proper ergonomics

Bariatric Vinyl Binder Lift Shown

Removable leg straps for superior versatility



Nylon

Vinyl

The Nylon Binder Lift received the JEMS 2013 Hot Product of the Year award and EMS World 2013 Innovation of the Year award. Constructed with Cordura Denier 1000 nylon and polypropylene web strapping for the handles. The removable leg straps are made from BioThane® and are 100% impervious to bodily fluids. This makes decontamination of the leg straps quick and easy. The core of the device is machine washable. Covered by a 5 year full warranty and has an expected use life of 10+ years.

BINDER LIFT™




The Vinyl Binder Lift received the JEMS 2014 Hot Product of the Year award and is our premier model. Constructed with Bondcote® vinyl coated polyester, and BioThane® for the handles and removable leg straps. Both materials used are 100% impervious to bodily fluid and are chemical resistant. To decontaminate simply wipe down with a sani-wipe and place back in the storage bag for next use. Covered by a 5 year full warranty and has an expected use life of 10+ years.

Available in two sizes:

STANDARD (chest range of 34-62 inches)

BARIATRIC (chest range of 58-82 inches)



Distributed by:
EA Medical, LLC • www.eamed.com
For more information or to schedule a demonstration contact:
Neal Zeid (314) 974-5915 email: nealrz@juno.com