

LUMALIER Disinfection Wall for PPE Doffing

Portable, UVC
Germicidal energy
wall for high-level
decontamination
and disinfection of
PPE equipment prior
to doffing.

A safe, advanced
engineering control
engineered to be
used as an adjunct to
existing doffing
protocols.

4-6 log pathogen
reduction on all
treated surfaces in
about one minute.



Safer Doffing = Increased Risk Reduction

The Lumalier PPE Disinfection Wall adds an important layer of infectious pathogen protection for your staff.

- There are no known pathogens that are resistant to a lethal dose of UVC germicidal energy...NONE.
- Portable, stackable, storable and transportable.
- Sets up / folds down in minutes.

Specifications:

- Total germicidal output ~**12,000 UVW / minute.**
~200 UVC watts per second @ one meter.
- The Lumalier PPE disinfection wall features twelve 55W Philips twin-tube non-proprietary twin-tube germicidal lamps, 44" total tube length per lamp, with a total of 528 inches of tube length.



- Green technology; lowest Hg = 5mg or less Mercury content.
- 9000 hours of continuous use. Philips twin-tube germicidal lamp.
- Simple 120V grounded 9' power cord.
- All aluminum and hospital-grade stainless-steel construction, with reflective back wall to direct and focus cleansing UV energy towards the targeted PPE equipment.
 - In closed, storage mode: the Lumalier PPE disinfection wall folds to a four-foot x six-foot x five-inch stackable unit.
 - In open, full disinfection mode: the Lumalier PPE disinfection wall stands six-feet tall x six-feet wide x three-feet deep.
- The Lumalier PPE disinfection wall folds quickly and is easily transportable. Wheels and handles help with transportation and set up.

Ballast: six non-proprietary Philips / Advance PowerVolt UV, self-regulating electronic ballast. 5year, non-prorated full replacement warranty.

Safety requirement:

UVC can cause temporary damage and discomfort to unprotected skin and eyes.

ALWAYS wear Personal Protective Equipment (PPE) for eyes and skin to prevent exposure (direct and / or indirect) to UVC germicidal energy.

**Lumalier Environmental Disinfection
Products By Evergreen UV LLC
Made in the USA**



Distributed by: EA Medical, LLC • 800-339-9393 • uvc@eamed.com
www.eamed.com • 10541 State Highway 81 • Canton, MO 63435

